



THIS LAYOUT IS REFLECTIVE OF THE LATEST SALES CONFIGURATION. ANY CHANGES TO THIS SALES CONFIGURATION MAY REQUIRE A REVISION TO THIS PROJECT PLAN.

SINKS, COUNTERTOPS AND ALL CASEWORK SHOWN IS SUGGESTED AND MUST BE DESIGNED SUPPLIED AND INSTALLED BY CUSTOMER/CONTRACTOR.

STRETCHER SIZE (6'-5" X 2'-7") SHOWN IS FOR REFERENCE ONLY. VERIFICATION AND COORDINATION BY CUSTOMER IS REQUIRED TO ENSURE PROPER TRANSPORT AND WORKFLOW ACCESS.



NOTE: THE UPS IS SUPPLIED AND DELIVERED TO CUSTOMER'S LOADING DOCK BY SIEMENS. CUSTOMER'S ELECTRICIAN IS RESPONSIBLE FOR MOVING FROM LOADING DOCK TO FINAL LOCATION AND COMPLETING ALL CONNECTIONS. THE UPS MUST NOT BE LOCATED IN A PATIENT ENVIRONMENT. SIEMENS PROJECT MANAGER WILL SCHEDULE UPS STARTUP PRIOR TO DELIVERY OF SIEMENS IMAGING EQUIPMENT.



Click to acquire Typical PDF technical data

SCALE: 1/8" = 1'-0"

PROJECT MANAGER: JOAS AGUILAR VALLEJO  
PHONE: (817) 366-5832  
EMAIL: JOAS.AGUILARVALLEJO@SIEMENS-HEALTHINEERS.COM

## JPS HEALTH NETWORK

1500 SOUTH MAIN STREET, FORT WORTH, TX 76104

**SIEMENS**

B 09/19/25 UPDATED SKETCH: CHANGES TO EQUIP RM / CABINETS PLACEMENT

SHEET OF 1 2

DATE: 09/19/25

























DRAWN BY: O. CARRILLO

REF: 30301324

IR LAB #4 2126 - ARTIS ICONO BIPLANE IR  
REFERENCE DOCUMENT - NOT FOR CONSTRUCTION

2411065

# EQUIPMENT LEGEND

NO	DESCRIPTION	SMS SYM	WEIGHT (LBS)	BTU/HR TO AIR	DIMENSIONS (INCHES)			REMARKS
					W	D	H	
①	CPC (CENTRAL POWER CONTROL)		3	N/A	9 5/8	4	2 1/2	ON CONTROL COUNTER
②	CONTROL INTERFACE BOARD		18	342	19 1/4	3 1/4	13 1/2	ON WALL UNDER COUNTER
③	KEYBOARD		2.2	342	17 1/2	6 1/8	2 1/8	MTD. UNDER COUNTER
④	24" CONTROL ROOM DISPLAY		17	191	21 3/4	8	20 1/2	ON COUNTER
⑤	INJECTOR WALL CONNECTION – CONTROL ROOM		18	342	20 1/4	4 3/4	13 3/8	ON WALL UNDER COUNTER
⑥	TABLE CONTROL MODULES		16	---	20	8 3/4	3 1/2	ON TABLE
⑦	DCS LARGE DISPLAY RAIL MOUNT W/ 1.12 METER EXTENSION RAILS		787	1,535	44	60 3/4	50 3/4	CEILING SUSPENDED
⑧	ARTIS ICONO BIPLANE FLOOR STAND W/ MOUNTING PLATE		1,477	682	---	---	---	C-ARM FLOOR MOUNTED
⑨	ARTIS ICONO BIPLANE CEILING STAND W/ LONGITUDINAL RAILS		2,443	682	---	---	---	C-ARM CEILING MOUNTED
⑩	PATIENT TABLE (MULTI-TILT)		1,200	682	---	---	---	TABLE FLOOR MOUNTED
⑪	UPPER BODY RADIATION SHIELD 4 M TRACK		196	---	---	---	---	TRACK MOUNTED
⑫	MAVIG LAMP		48	---	---	---	---	
⑬	GENERATOR 1 (ACX)		551	4,095	23 5/8	23 1/4	63 1/2	FLOOR MOUNTED
⑭	GENERATOR 2 (ACX)		551	4,095	23 5/8	23 1/4	63 1/2	FLOOR MOUNTED
⑮	CABLE CABINET		265	---	31 1/2	17 1/8	87	FLOOR MOUNTED
⑯	SYSTEM CONTROL CABINET		970	13,649	39 1/2	25 1/2	74 3/4	FLOOR MOUNTED
⑰	ATIS IMAGE SYSTEM		772	6,483	39 1/2	25 1/2	74 3/4	FLOOR MOUNTED
⑱	TUBE COOLING UNIT (PLANE A)		80	15,355	16 1/2	28 1/4	19 1/4	FLOOR MOUNTED
⑲	TUBE COOLING UNIT (PLANE B)		80	15,355	16 1/2	28 1/4	19 1/4	BRACKET MOUNTED
⑳	EATON 9355 15KVA UPS AND BATTERY		755	8,134	12 3/4	33 1/2	47 3/4	SEE MFG REQUIREMENTS
㉑	EATON 9355 OUTPUT TRANSFORMER CABINET		490	---	20	34 1/8	66	SEE MFG REQUIREMENTS
㉒	EATON 9355 REMOTE MONITORING DEVICE		0.5	---	6	1	3	SEE MFG REQUIREMENTS
㉓	VITALINQ INTERCOM CONSOLE		---	---	10	10	5	INSTALLED BY CUSTOMER
㉔	VITALINQ MICROPHONE		---	---	---	---	---	INSTALLED BY CUSTOMER

RESTRICTIONS:  
THIS LAYOUT MEETS SYSTEM REQUIREMENTS WITH NO RESTRICTIONS.

## Authorization:

I have reviewed the below referenced plan(s) and understand and agree with the proposed room configuration and equipment arrangement. By signing below, I agree that I have the authority to sign on behalf of the customer. I have been fully informed by Siemens about the limitations (if applicable) that exist regarding the proposed equipment installation. Siemens has explained that the listed limitations may interfere with certain functions of the equipment. We agree to assume responsibility to inform employees, agents, contractors, users or any other persons handling this equipment of these limitations and will assure that the equipment is handled in compliance with these parameters.

Signature: \_\_\_\_\_ Date: \_\_\_\_\_

Name Printed: \_\_\_\_\_ Title: \_\_\_\_\_

PROJECT MANAGER: JOAS AGUILAR VALLEJO  
PHONE: (817) 366-5832  
EMAIL: JOAS.AGUILARVALLEJO@SIEMENS-HEALTHINEERS.COM

## JPS HEALTH NETWORK

1500 SOUTH MAIN STREET, FORT WORTH, TX 76104

B 09/19/25 UPDATED SKETCH: CHANGES TO  
EQUIP RM / CABINETS PLACEMENT

SHEET OF  
2 2

DATE:  
09/19/25

DRAWN BY:  
O. CARRILLO

REF:  
30301324

IR LAB #4 2126 - ARTIS ICONO BIPLANE IR  
REFERENCE DOCUMENT - NOT FOR CONSTRUCTION

# SIEMENS

2411065

## ARCHITECTURAL NOTES

- 1) ALL PRELIMINARY EQUIPMENT LAYOUTS SUBMITTED BY SIEMENS HEALTHCARE ARE BASED ON THE RECOMMENDED SPACE NECESSARY FOR THE OPERATION AND SERVICEABILITY OF THE EQUIPMENT BEING PROPOSED. SIEMENS WILL NOT SUBMIT AN EQUIPMENT LAYOUT THAT IS NOT IN THE BEST INTEREST OF BOTH THE CUSTOMER AND SIEMENS. ALL EQUIPMENT LAYOUTS ARE BASED EITHER ON AN ACTUAL SITE SURVEY OR ARCHITECTURAL DRAWINGS SUPPLIED TO SIEMENS. SIEMENS WILL NOT BE RESPONSIBLE FOR ANY ALTERATIONS THAT ENCRANCH WITHIN DESIGNATED SAFETY AND SERVICE CLEARANCE ZONES AS INDICATED ON DRAWINGS (I.E., PIPE CHASES, VENTILATION DUCTS, CASEWORK, AND SOFFITS, ETC.) MADE BY THE CUSTOMER OR REQUIRED BY A CUSTOMER'S ARCHITECTURAL FIRM ONCE PRELIMINARY DRAWINGS HAVE BEEN SUBMITTED AND APPROVED. DO NOT ALTER ANY SPECIFICATIONS AND/OR DIMENSIONS WITHOUT CONTACTING AND RECEIVING WRITTEN CONFIRMATION FROM SIEMENS PROJECT MANAGER.
- 2) SIEMENS HEALTHCARE IS NOT AN ARCHITECTURAL OR ENGINEERING FIRM. DRAWINGS SUPPLIED BY SIEMENS ARE NOT CONSTRUCTION DRAWINGS. THEREFORE, THESE DRAWINGS ARE TO BE USED ONLY FOR INFORMATION TO COMPLEMENT ACTUAL CONSTRUCTION DRAWINGS AVAILABLE FROM A CUSTOMER APPOINTED ARCHITECTURAL REPRESENTATIVE OR A CUSTOMER'S ENGINEERING DESIGN GROUP. THE CUSTOMER'S ARCHITECT AND GENERAL CONTRACTOR SHALL BE ULTIMATELY RESPONSIBLE FOR COMPLIANCE WITH ALL APPLICABLE CODES AND PROFESSIONAL DESIGN REQUIREMENTS INCLUDING OSHA/NEC SAFETY CLEARANCE REQUIREMENTS IN ADDITION TO SIEMENS-REQUIRED SAFETY/SERVICE CLEARANCES SHOWN.
- 3) THE CUSTOMER IS RESPONSIBLE FOR ALL ROOM AND AREA PREPARATION COSTS, PROFESSIONAL FEES, PERMITS, REPORTS, AND INSPECTION FEES.
- 4) EQUIPMENT WARRANTIES, EXPRESSED OR IMPLIED ON THE PART OF SIEMENS SHALL BE CONTINGENT UPON STRICT COMPLIANCE WITH THE ARCHITECTURAL, STRUCTURAL, ELECTRICAL, MECHANICAL AND RECOMMENDATIONS AND REQUIREMENTS CONTAINED IN THESE DRAWINGS, UNLESS SPECIFIED OTHERWISE.
- 5) ALL DIMENSIONS SHOWN ARE FROM FINISHED SURFACES UNLESS SPECIFIED OTHERWISE.
- 6) THIS DRAWING DOES NOT PROVIDE RADIATION SHIELDING REQUIREMENTS FOR X-RAY AND ASSOCIATED EQUIPMENT. THE CUSTOMER IS RESPONSIBLE FOR CONSULTING WITH A REGISTERED RADIATION PHYSICIST. ACTUAL PROTECTION REQUIREMENTS SHALL BE SPECIFIED BY A REGISTERED RADIATION PHYSICIST AT CUSTOMER'S ENGAGEMENT AND EXPENSE. RESPONSIBILITY FOR ALL INFORMATION AS TO THE ROOM LOCATION, USE, AND NUMBER OF ANTICIPATED EXAMINATIONS TO BE PERFORMED PER TIME PERIOD SHALL BE PROVIDED TO THE PHYSICIST BY THE CUSTOMER. THE CUSTOMER SHALL FURTHER TAKE ALL RESPONSIBILITY IN THE COMMUNICATION AND COORDINATION OF ACTIVITIES OF THE RADIATION PHYSICIST AND THE ARCHITECTURAL REPRESENTATIVE.
- 7) SIEMENS HEALTHCARE SHALL BE RESPONSIBLE FOR SIEMENS EQUIPMENT INSTALLATION, CALIBRATION, CONNECTION AND INSTALLATION OF SIEMENS PROVIDED CABLES. THE CUSTOMER/ELECTRICAL CONTRACTOR IS RESPONSIBLE FOR TERMINATIONS OF CUSTOMER/ELECTRICAL CONTRACTOR-SUPPLIED CABLES TO SIEMENS EQUIPMENT. IN THE EVENT THAT SPECIFIC TRADE RULES OR LICENSE REQUIREMENTS PROHIBIT THIS, THE CUSTOMER SHALL INITIATE THE SERVICES OF APPROVED OTHER CONTRACTORS AND PAY FOR SELECTED, APPROVED PARTIES TO PERFORM THIS WORK WITH SUPERVISION PROVIDED BY SIEMENS. CALIBRATION WHEN ACCOMPLISHED OUTSIDE OF NORMAL INSTALLATION SEQUENCES DUE TO CONTRACTOR OR TRADE RULE ACTIONS OR REQUIREMENTS SHALL BE SUPPORTED BY, CHARGED TO, AND ACCEPTED BY THE CUSTOMER AS AN ADDITIONAL INSTALLATION EXPENSE.
- 8) THE CUSTOMER SHALL COORDINATE WITH SIEMENS PROJECT MANAGER THE LOCATIONS AND TRAVEL OF ALL ANCILLARY EQUIPMENT TO BE CEILING OR WALL MOUNTED (I.E.: O.R. LIGHTS, MEDICAL GAS COLUMNS, PHYSIOLOGICAL MONITORING INJECTORS, CRT PLATFORMS, SPRINKLER HEADS, SMOKE DETECTORS, ELECTRICAL OUTLETS, HVAC GRILLES, SPEAKERS, AND GENERAL ROOM LIGHTING, ETC.).
- 9) THE GENERAL CONTRACTOR/CUSTOMER SHALL BE RESPONSIBLE FOR ALL FINAL PAINT, TOUCH-UP AND ANY COSMETIC OR TRIM WORK WHICH NEEDS TO BE OR IS REQUIRED TO BE COMPLETED AFTER THE INSTALLATION OF THE SIEMENS EQUIPMENT AND ANY ASSOCIATED SUPPORT APPARATUS.
- 10) CUSTOMER/CONTRACTOR MUST ASSIST SIEMENS INSTALLERS WITH INSTALLATION OF EQUIPMENT ABOVE 14'-0". REFER TO THE ELECTRICAL NOTES ON SIEMENS SHEET E-101 FOR MORE DETAILS.

FOR ADDITIONAL PLANNING RELATED QUESTIONS,  
PLEASE REFER TO TYPICAL DRAWING SET OR  
SIEMENS PROJECT MANAGEMENT.